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IHEJ-SINR RC RM-J [Radiology: SINR]

Report Submission to Report Manager

Report Creators send SR objects to Report Manager. Report actors should attempt the following cases. Some may not apply because the actors may not support that feature.

1. Simple text report with no image references.
2. Simple text report with single image reference.
3. Simple text report with several image references
4. Numeric report

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Test script :

Step	Trans.	Msg Type	Return	From Actor	To Actor	Step Description	Opt.
1	NULL	C-Store		RC	RM	Report Creator sends report to Report Manager	R

IHEJ-SINR RC RRD-J [Radiology: SINR]

Tests of Creator/Reader

Report Creators send SR objects to Report Readers (possibly through a Report Repository). Report Reader and renders report. Images that are referenced should be retrieved and rendered. Report actors should attempt the following cases. Some may not apply because the actors may not support that feature.

1. Simple text report with no image references.
2. Simple text report with single image reference.
3. Simple text report with several image references
4. Numeric report
5. Report amended by Report Creator

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Test script :

Step	Trans.	Msg Type	Return	From Actor	To Actor	Step Description	Opt.
1	NULL			RC	RRD	Report Creator sends report to Report Reader (not Manager)	R
2	NULL			RRD	--	Report Reader renders report (+ images)	R

IHEJ-SINR RM RRD-J [Radiology: SINR]

Tests of Report Manager/Reader

Report Manager modifies and sends SR objects to Report Readers (possibly through a Report Repository). Report Reader and renders report. Images that are referenced should be retrieved and rendered. Report actors should attempt the following cases. Some may not apply because the actors may not support that feature.

1. Simple text report with no image references.
2. Simple text report with single image reference.
3. Simple text report with several image references
4. Numeric report

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Test script :

Step	Trans.	Msg Type	Return	From Actor	To Actor	Step Description	Opt.
1	NULL	C-store		RM	RRD	Report Manager sends report to Report Reader (not Manager)	R
2	NULL			RRD	--	Report Reader renders report (+ images)	R

IHEJ-SINR RM RRP-J [Radiology: SINR]

Test the submission of reports to report repository.

This test is mainly used for testing the storage between the report manager and the report repository actor.

This test makes sense if the report manager and the report repository are not grouped on the same system.

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Test script :

Step	Trans.	Msg Type	Return	From Actor	To Actor	Step Description	Opt.
11	RAD-25	C-Store		RM	RM	Change the SR Status to verified	R
12	RAD-25	C-Store		RM	RRP	store the SR	R

IHEJ-SINR RRD RRP-J [Radiology: ARI, SINR]

ARI SR Objects

Report Reader should query Report Repository for SR objects. Report Reader should demonstrate query keys as listed in TF.

SR objects should be retrieved and rendered. If reports contain references to images, images should be retrieved and displayed (if supported).

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Test script :

Step	Trans.	Msg Type	Return	From Actor	To Actor	Step Description	Opt.
100	RAD-26	Query SR Objects(C-FIND)	Responses (C-FIND)	RRD	RRP	Query Report Repository forSR objects	R
210	RAD-27	Retrieve SR Objects(C-STORE)	Store SR Objects(C-STORE)	RRD	RRP	Retrieve SR Objects	R

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IHEJ-CT TS TC [IT Infrastructure: CT, (Radiology: SEC)]

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Integration Profiles / Actors Concerned by the test :

<u>Integration Profile</u>	<u>Actor</u>	<u>Profile Option</u>	<u>Option</u>	<u>Meta Test</u>
CT	TIME_CLIENT	NONE	R	
CT	TS	NONE	R	

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"IHEJ-LSWF test" was divided like the following list.

1. IHEJ-LAB_ADT (include a part of "LIR")
2. IHEJ-LAB_ORDER
3. IHEJ-LAB_OF_ORT
4. IHEJ-LAB_FILLER_ORDER_MANAGEMENT(not to do)
5. IHEJ-LAB_OF_AM
6. IHEJ-LAB_AM_OF

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7. IHEJ-LAB_AM_LD
8. IHEJ-LAB_LD_AM

IHEJ-LAB ADT [Laboratory]

Laboratory Scheduled Workflow Testing the Patient Registration

ADT regist a patient information to OP,OF,ORT

This test include a part of "LIR"(Change Patient Name)

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Integration Profiles / Actors Concerned by the test :

<u>Integration Profile</u>	<u>Actor</u>	<u>Profile Option</u>	<u>Option</u>	<u>Meta Test</u>
LSWF	ADT	NONE	R	
LSWF	OF	NONE	R	
LSWF	OP	NONE	R	
LSWF	ORT	NONE	O	

Test script :

Step	Trans.	Msg Type	Return	From Actor	To Actor	Step Description	Opt.
1	RAD-12	ADT^A08	ACK^A08	ADT	OP	New Patient Registration	R
2	RAD-12	ADT^A08	ACK^A08	ADT	OF	New Patient Registration	R
4	RAD-12	ADT^A08	ACK^A08	ADT	OP	Change Patient Name (a part of LIR)	O
5	RAD-12	ADT^A08	ACK^A08	ADT	OF	Change Patient Name (a part of LIR)	O

IHEJ-LAB AM LD [Laboratory]

Laboratory Work Order Step

Automation Manager place a new oeder to the Laboratory Device.

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Integration Profiles / Actors Concerned by the test :

<u>Integration Profile</u>	<u>Actor</u>	<u>Profile Option</u>	<u>Option</u>	<u>Meta Test</u>
LDA	AM	NONE	R	
LDA	ANALYZER	NONE	R	
LDA	PRE_POST_PROCESSOR	NONE	R	

Test script :

Step	Trans.	Msg Type	Return	From Actor	To Actor	Step Description	Opt.
1	LAB-21	OML^O33	ORL^O34	AM	ANALYZER	New Order	R
2	LAB-21	OML^O33	ORL^O34	AM	PRE_POST_PROCESSOR	New Order	R

IHEJ-LAB AM OF [Laboratory]

Laboratory Test Results Management

This test corresponds to transaction LAB-5 of IHE Laboratory Technical Framework. The actors using this transaction are the Order Filler and the Automation Manager.

The Order Filler manages the test results notified by the Automation Manager.

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Integration Profiles / Actors Concerned by the test :

<u>Integration Profile</u>	<u>Actor</u>	<u>Profile Option</u>	<u>Option</u>	<u>Meta Test</u>
LSWF	AM	NONE	R	
LSWF	OF	NONE	R	

Test script :

Step	Trans.	Msg Type	Return	From Actor	To Actor	Step Description	Opt.
1	LAB-5	OUL^R22,R23,R24	ACK^R22,R23,R24	AM	OF	Results are available	R

IHEJ-LAB-FILLER ORDER MANAGEMENT [Laboratory]

Laboratory Scheduled Workflow Test Filler Order Management

This test doesn't used in the "IHE-J 2006 Connectathon".

It was reserved for future.

The purpose here is to test the transaction LAB-2

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Integration Profiles / Actors Concerned by the test :

<u>Integration Profile</u>	<u>Actor</u>	<u>Profile Option</u>	<u>Option</u>	<u>Meta Test</u>
LSWF	OF	NONE	R	
LSWF	OP	NONE	R	

Test script :

Step	Trans.	Msg Type	Return	From Actor	To Actor	Step Description	Opt.
1	LAB-2	OML^O33		OF	OP	OML Send order/service number: ORC-1 = SN	O
2	LAB-2	ORL^O34		OP	OF	ORL acknowledgement Number assigned: ORC-1 = NA Unable to accept order: ORC-1 = UA	O

IHEJ-LAB LD AM [Laboratory]

Laboratory Work Order Step

Laboratory Device send a query to the Automation Manager(with barcode information).

Automation Manager response to the Laboratory Device with order.

Laboratory Device send observation result or specimen status change to the Automation Manager.

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Integration Profiles / Actors Concerned by the test :

<u>Integration Profile</u>	<u>Actor</u>	<u>Profile Option</u>	<u>Option</u>	<u>Meta Test</u>
LSWF	AM	NONE	R	
LSWF	OF	NONE	R	

Test script :

Step	Trans.	Msg Type	Return	From Actor	To Actor	Step Description	Opt.
1	LAB-4	OML^O33	ORL^O34	OF	AM	New order entered at OF	R
2	LAB-4	OML^O33	ORL^O34	OF	AM	Order Modify	O
3	LAB-4	OML^O33	ORL^O34	OF	AM	Cancel Order	O

IHEJ-LAB OF ORT [Laboratory]

Order Filler to ORT(OP) transaction

Oder Filler provides notification to the Order Result Tracker (or Order Placer) for specimen arrival, observation results.

This test consists 3 transactions.

1. Specimen Arrival
2. Observation Results (Result status A or R)
3. Observation Results (Result status F)

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Integration Profiles / Actors Concerned by the test :

<u>Integration Profile</u>	<u>Actor</u>	<u>Profile Option</u>	<u>Option</u>	<u>Meta Test</u>
LSWF	OF	NONE	R	
LSWF	ORT	NONE	R	

Test script :

Step	Trans.	Msg Type	Return	From Actor	To Actor	Step Description	Opt.
1	LAB-3	OUL^R22	ACK^R22	OF	ORT	After reception and acceptance of a specimen, inform the ORT(OP)	R
2	LAB-3	OUL^R22	ACK^R22	OF	ORT	Acquisition of technically validated result (Result status A or R)	O
3	LAB-3	OUL^R22	ACK^R22	OF	ORT	Clinical validation of results (Change Status F)	R
4	LAB-3	OUL^R22	ACK^R22	OF	ORT	Add & delete some test items	O

IHEJ-LAB ORDER [Laboratory]

Laboratory Scheduled Workflow Testing the Order

Oder Placer place a new oeder to the Oeder Filler.

This test consists 3 transactions.

1. New Order without barcode information
2. Cancel Order
3. New Order with barcode information

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Integration Profiles / Actors Concerned by the test :

<u>Integration Profile</u>	<u>Actor</u>	<u>Profile Option</u>	<u>Option</u>	<u>Meta Test</u>
LSWF	OF	NONE	R	
LSWF	OP	NONE	R	

Test script :

Step	Trans.	Msg Type	Return	From Actor	To Actor	Step Description	Opt.
1	LAB-1	OML^O33		OP	OF	New Order without Barcode information	O
2	LAB-1	ORL^O34		OF	OP	Response to the order	O
3	LAB-1	OML^O33		OP	OF	Cancel Order	O
4	LAB-1	ORL^O34		OF	OP	Response to the order	O
5	LAB-1	OML^O33		OP	OF	New Order with Barcode information	R
6	LAB-1	ORL^O34		OF	OP	Response to the order	R
7	LAB-1	OML^O33		OP	OF	Add & Delete some test items	R
8	LAB-1	ORL^O34		OF	OP	Response to the order	R

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C-ID CATH [Cardiology: CATH]

Image Display Actor Displays CATH Images

The purpose of this test is to verify the Image Display of variety of CATH image sets from various vendors.

Integration Profiles / Actors Concerned by the test :

<u>Integration Profile</u>	<u>Actor</u>	<u>Profile Option</u>	<u>Option</u>	<u>Meta Test</u>
CATH	ID	NONE	R	
CATH	IM	NONE	R	
CATH	MOD	NONE	R	

Test script :

Step	Trans.	Msg Type	Return	From Actor	To Actor	Step Description	Opt.
20	RAD-8	C-STORE (Images Stored)		MOD	IM	Send CATH images to Central Archive	R
40	RAD-14	Query Images (C-FIND)	Response	ID	IM	Query the IM for all XA images.	R
60	CARD-4	Retrieve Images (C-MOVE)	Store Images (C-STORE)	ID	IM	Retrieve all available images.	R
80	NULL	Display Images	verify the image encoding and patient demographics are accurate.	ID	ID	Display every image loop of every study	R

C-ID ECHO [Cardiology: ECHO]

Image Display Actor Displays ECHO Images

The purpose of this test is to verify the Image Display of variety of ECHO image sets from various vendors.

The Image Display shall display every image loop of every study and verify the image encoding and patient demographics are accurate.

We are not testing the query and retrieve of images, we are testing the ability of the ID to correctly render the images. The central archive can be used to store set of images from different vendors.

Integration Profiles / Actors Concerned by the test :

<u>Integration Profile</u>	<u>Actor</u>	<u>Profile Option</u>	<u>Option</u>	<u>Meta Test</u>
ECHO	ID	NONE	R	
ECHO	IM	NONE	R	
ECHO	MOD	NONE	R	

Test script :

Step	Trans.	Msg Type	Return	From Actor	To Actor	Step Description	Opt.
11	NULL	C-STORE		MOD	IM	Send ECHO images to Central Archive	O
40	RAD-14	Query Images (C-FIND)	Responses (C-FIND)	ID	IM	Query the Central Archive for all US images.	R
60	CARD-4	Retrieve Images (C-MOVE)	Store Images (C-STORE)	ID	IM	Retrieve all available images.	R
80	NULL	Render retrieved images		ID	ID	Display Images	R

C-ID ECHO STRESS [Cardiology: ECHO]

Image Display Actor Displays Stress ECHO Images from Modality

The purpose of this test is to verify the Image Display of variety of ECHO STRESS image sets from various vendors. The only required actor for this test is the IMAGE DISPLAY.

It is recommended that this test is combined with C-MOD-STRESS to reduce the number of tests which must be run.

Integration Profiles / Actors Concerned by the test :

<u>Integration Profile</u>	<u>Actor</u>	<u>Profile Option</u>	<u>Option</u>	<u>Meta Test</u>
ECHO	ID	STRESS_ECHO	R	
ECHO	IM	NONE	O	

Test script :

Step	Trans.	Msg Type	Return	From Actor	To Actor	Step Description	Opt.
40	RAD-14	Query Images (C-FIND)	Responses (C-FIND)	ID	IM	Query the Central Archive for US stress echo images.	R
60	CARD-4	Retrieve Images (C-MOVE)	Store Images (C-STORE)	ID	IM	Retrieve all available images.	R
80	NULL	Render retrieved images		ID	ID	Display stress echo Images	R

C-IM CATH [Cardiology: CATH]

Image Manager Accepts/Stores CATH Images

The purpose of this test is to verify that an Image Manager can accept and store all real world vendor images which are available.

Integration Profiles / Actors Concerned by the test :

<u>Integration Profile</u>	<u>Actor</u>	<u>Profile Option</u>	<u>Option</u>	<u>Meta Test</u>
CATH	IM	NONE	R	

Test script :

Step	Trans.	Msg Type	Return	From Actor	To Actor	Step Description	Opt.
20	CARD-2			--	IM	Either query for or have the Connectathon Image Repository send all available XA images to the Image Manager.	R
40	NULL			--	--	Verify that all images are received and can be queried and retrieved from the IM under test. Verify that the objects are stored as full-fidelity objects (ie., no DICOM data elements are deleted).	R

C-IM ECHO [Cardiology: ECHO]

Image Manager Accepts/Stores ECHO Images

The purpose of this test is to verify that an Image Manager can accept and store all real world vendor images which are available.

Integration Profiles / Actors Concerned by the test :

<u>Integration Profile</u>	<u>Actor</u>	<u>Profile Option</u>	<u>Option</u>	<u>Meta Test</u>
ECHO	IM	NONE	R	

Test script :

Step	Trans.	Msg Type	Return	From Actor	To Actor	Step Description	Opt.
20	CARD-2			--	IM	Either query for or have the Connectathon Image Repository send all available US images to the Image Manager.	R
40	NULL			--	--	Verify that all images are received and can be queried and retrieved from the IM under test. Verify that the objects are stored as full-fidelity objects (ie., no DICOM data elements are deleted).	R

C-IM SCHED ECHO [Cardiology: ECHO]

ECHO Intermittently Connected Devices

The purpose of this test is to verify that the Storage Commitments are handled correctly with Intermittently connected modalities, especially Echo (see CARD-3).

Specifically, this is Use Case E2. This is the simplest case where the Patient has already been Registered, the Order has already been created, and the Procedure was already been scheduled.

Integration Profiles / Actors Concerned by the test :

<u>Integration Profile</u>	<u>Actor</u>	<u>Profile Option</u>	<u>Option</u>	<u>Meta Test</u>
ECHO	IM	NONE	R	
ECHO	MOD	NONE	R	

Test script :

Step	Trans.	Msg Type	Return	From Actor	To Actor	Step Description	Opt.
20	NULL			--	--	Admit a patient called IMECHO^MOD scompany name, for example, IMECHO^PHILIPS.	R
40	NULL			--	--	Verify that all images are received and can be queried and retrieved from the IM under test. Verify that the objects are stored as full-fidelity objects (ie., no DICOM data elements are deleted).	R
60	NULL			--	--	At a DSS/OF schedule this procedure.	R
80	RAD-4			OF	IM	DSS/OF notifies IM of produce. (Procedure Scheduled)	R
100	RAD-5			MOD	OF	Modality Query for Echo Worklist.	R
120	NULL			--	--	Literally disconnect the Echo Modality from the network	R
160	NULL			--	--	Reconnect the Echo system to the network	R

180	NULL			--	--	Perhaps consider trying to complete CARD-1, CARD-2, and RAD-7 and then, with excellent timing, disconnecting the echo system from the network prior to CARD-3. Attempt to force the IM Storage Commit response to fail, then retry upon the next association from the Echo Modality.	R
200	NULL			--	--	Note that the following 4 transactions do not have to occur in exact order.	R
220	CARD-1			MOD	PPSM	MPPS In Progresss from Echo Modality.	R
240	CARD-2			MOD	ID	From Echo modality send study to the IM. (C-store)	R
260	RAD-7			MOD	PPSM	MPPS Complete from Echo Modality	R
300	NULL			--	--	Verify that all images within the study are sent.	R
320	NULL			--	--	Verify that the images created at the modality contain the OF Generated Study Instance UID.	R

C-IM SCHED ECHO NEW [Cardiology: ECHO]

Echo Use Case E5

The purpose of this test is to verify that the Storage Commitments are handled correctly with Intermittently connected modalities, especially Echo (see CARD-3).

Specifically, this is Use Case E5. This is the **worst case** where the Patient has not been Registered, the Order was not created, and the Procedure was not scheduled.

Integration Profiles / Actors Concerned by the test :

<u>Integration Profile</u>	<u>Actor</u>	<u>Profile Option</u>	<u>Option</u>	<u>Meta Test</u>
ECHO	IM	NONE	R	
ECHO	MOD	NONE	R	

Test script :

Step	Trans.	Msg Type	Return	From Actor	To Actor	Step Description	Opt.
20	NULL			MOD	IM	Literally disconnect the Echo Modality from the network.	R
60	NULL			--	--	Reconnect the Echo system to the network	R
80	NULL			--	--	Note that the following 4 transactions do not have to occur in exact order	R
100	CARD-1			MOD	PSPM	MPPS In Progresss from Echo Modality	R
120	CARD-2			MOD	IM	From Echo modality send study to the IM. (C-store)	R
140	RAD-7			MOD	PSPM	MPPS Complete from Echo Modality	R
180	NULL			--	--	Verify that all images within the study are sent	R
200	NULL			--	--	Verify that the images created at the modality contain the MODALITY Generated Study Instance UID.	R

C-IMOF [Cardiology: CATH, ECHO]

Image Manager Updates DICOM Objects Whenever Patient Update is Received

The purpose of this test is to verify that the IM updates ALL image objects even if a Patient Update message is received in the middle of a procedure. That is, images which have already been received and images which have yet to be received must all be updated accurately.

It is recommended that IM actors find another IM partner and swap this test role with them, ie., run this test twice successively alternating roles. Note that this test only needs to be executed once between the CATH and ECHO Integration Profiles.

Integration Profiles / Actors Concerned by the test :

<u>Integration Profile</u>	<u>Actor</u>	<u>Profile Option</u>	<u>Option</u>	<u>Meta Test</u>
CATH	IM	NONE	R	
ECHO	IM	NONE	R	

Test script :

Step	Trans.	Msg Type	Return	From Actor	To Actor	Step Description	Opt.
20	NULL			OF	OF	Admit a patient called IMUPDATE^IMu company name, for example, IMUPDATE^EMAGEON.	R
21	NULL			OF	OF	On some other IM or MESA test tool, create an image set of the same patient name and ID.	R
22	CARD-2			IM	IM	From some device, modality or other image sending option, begin to send a large study to the IM. (C-store)	R
23	RAD-12			OF	IM	From the OF (or MESA RIS Mail) send a Patient Update message to the IM, in the middle of sending the image set, if at all possible. Change the patient name to, eg., IMUPDATE^EMAGEON2 .	O
24	NULL			IM	IM	Verify that all images within the study are sent.	R
25	RAD-14			ID	IM	From some other device query the IM under test for the original patient Name. There should be no images found.	R
26	RAD-14			ID	IM	From some other device query the IM under test for the update patient Name (with 2). Verify that all images are found.	R

C-MOD-ECHO STRESS

Verify Encoding of ECHO Stress Images

The purpose of this test is to verify the Modality of ECHO STRESS image sets are properly encoding the required STRESS ECHO DATA ELEMENTS. The only required actor for this test is the ECHO MODALITY.

It is recommended that this test is combined with C-ID-STRESS to reduce the number of tests which must be run.

Note that it is recommended that the modality arrives at the Connectathon with these images PRE-OBTAINED and PRE-LOADED or you will waste a great deal of time!

Integration Profiles / Actors Concerned by the test :

<u>Integration Profile</u>	<u>Actor</u>	<u>Profile Option</u>	<u>Option</u>	<u>Meta Test</u>
ECHO	MOD	NONE	R	

Test script :

Step	Trans.	Msg Type	Return	From Actor	To Actor	Step Description	Opt.
20	NULL			MOD	--	Configure the Echo Modalities to send to an appropriate receiving device.	R
40	CARD-2			MOD	IM	Send a set of STRESS ECHO images which demonstrate as many as possible of the enumerated values of the following elements: Stage Number (0008,2122), View Number (0008,2128), Stage Code Sequence (0040,000A), and View Code Sequence (0008,2240) (see IHE Cardiology Technical Framework Volume 2 section 4.2.3).	R
60	NULL			--	ID	Either display every image loop of every study and verify the image encoding and patient demographics are accurate or manually edit the image headers to verify.	RO

C-MOD DSS/OF C3 C5 [Cardiology: CATH]

CATH Use Cases C3 and C5 for Modality and DSS/OF

Ideally, there would be a hemodynamic (HD) and an angio (XA) modality available for these tests, however, that will most likely not be the case. If there is no HD modality available, execute with only a single modality for a cath lab. In the real world, this is analogous to the hemo system being non-IHE compliant.

This test primarily tests the DSS/OF ability to handle the complexities of CARD-1 transaction. Specifically, this test tests CATH Use Cases C3 and C5 ability to create Requested Procedures and SPSs u201con the flyu201d with the correct Study Instance UID from the Modality.

Modality performs an unscheduled procedure. Order Filler dynamically creates a Requested Procedure/Scheduled Procedure Step for that modality and any other modalities in that lab. Image Manager may or may not also participate. The MPPS Mgr can be combined with either the DSS/OF or the IM in these test cases, but only one MPPS Mgr should be active at one time and it should be recorded which MPPS Mgr was used.

Integration Profiles / Actors Concerned by the test :

<u>Integration Profile</u>	<u>Actor</u>	<u>Profile Option</u>	<u>Option</u>	<u>Meta Test</u>
CATH	MOD	NONE	R	
CATH	OF	NONE	R	

Test script :

Step	Trans.	Msg Type	Return	From Actor	To Actor	Step Description	Opt.
20	NULL			--	--	In the descriptions to follow, SYS_NAME is a combination of the Modality(s) and DSS/Order Filer systems under test	R
40	NULL			--	--	For example: GEXAPHILIPSHDKODAKOF	R
60	NULL			--	--	Find and ADT and Order Placer to help with this test. This may be the MESA RIS Mail	R
80	NULL			--	--	Use the ADT system to register a single inpatient: SYS_NAME^JK	R
100	NULL			--	--	This begins the Un-ordered/unscheduled case	R
120	CARD-1			MOD	PPSM	HD or XA MPPS In Progress: SYS_NAME^JK	R
140	NULL			--	--	Verify on the DSS/OF that a RP/SPS has been created for this procedure using the modality provided Study Instance UID.	R

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C-MOD IM CATH [Cardiology: CATH]

CATH Profile: Modality/IM Communication

C-MOD_IM_CATH tests the interactions between Acquisition Modalities and Image Managers. The specific modalities included in this test are XA, HD or IVUS.

Modality performs different procedures (with different MPPS cases) and communicates those results with the Image Manager. This test assumes that the MPPS Mgr is associated with the IM. The IM/MMPS Mgr does not need to forward the MPPS messages to the OF for this test.

If you find a DSS/OF that wishes to participate, run the C-MOD_IM_DSS/OF test instead of this one. Note that this test case does NOT test the complexities of the DSS/OF for the CARD-1 transaction. See the C-MOD_IM_DSS/OF tests.

For evaluation, the Project Manger should visit the Image Manager and ask for evidence that all imaging steps were performed as requested. The Image Manager should so logs or user interface functions that indicate MPPS messages were received (IN PROGRESS, COMPLETE, DISCONTINUED). Image Managers with embedded displays should display images from the modality. Image Managers without displays can show evidence the images were stored (database entries, other GUI) or can direct the Project Manager to an external (other) Image Display.

Integration Profiles / Actors Concerned by the test :

Integration Profile	Actor	Profile Option	Option	Meta Test
CATH	IM	NONE	R	
CATH	MOD	NONE	R	

Test script :

Step	Trans.	Msg Type	Return	From Actor	To Actor	Step Description	Opt.
20	NULL			--	--	In the descriptions to follow, SYS_NAME is a combination of the Modality and Image Mgr systems under test	R
30	NULL			--	--	For example: PHILPISXASIEMENSIM	R
40	NULL			--	--	Find an ADT, Order Placer, DSS trio to help with this test. This may be the MESA RIS Mall.	R
50	NULL			--	--	Use the ADT system to register 5 inpatients: SYS_NAME^QR, SYS_NAME^ST, SYS_NAME^UV, SYS_NAME^WX, SYS_NAME^YZ	R
60	NULL			--	--	Use the Order Placer and place single orders for QR, ST, UV, WX	R
70	NULL			--	--	Use the Order Placer and place two orders for YZ that can be used for the Group Case	R

80	NULL			--	--	Find a OF/MWL provider and schedule 1 SPS or SYS_NAME^QR, ^ST, ^UV, ^WX,	R
90	NULL			--	--	Find a OF/MWL provider and schedule 2 SPSs for SYS_NAME^YZ	R
240	NULL			--	--	Step 3: This begins the PPS Abandoned Case for patient SYS_NAME^UV; select this SPS and abandon the step. Abandon by sending an MPPS In Progress, then MPPS Discontinued.	R
260	CARD-1			MOD	PPSM	MPPS for Abandoned Case: SYS_NAME^UV	R
280	NULL			--	--	Step 4: This begins the PPS Append Case for patient SYS_NAME^WX; perform one (scheduled) PPS for SYS_NAME^WX, and then append the same type of PPS	R
300	CARD-1			MOD	PPSM	MPPS for first step of Append Case: SYS_NAME^WX	R
320	CARD-2			MOD	ID	C-Store	R
340	CARD-1			MOD	PPSM	MPPS for Append PPS: SYS_NAME^WX	R
360	CARD-2			MOD	IM	C-Store	R
380	NULL			--	--	Step 5: This begins the PPS Group Case for patient SYS_NAME^YZ; modality groups two SPS in to one PPS. Note that this step is optional for CATH, only complete if the modality is capable of grouping procedures.	O
400	CARD-1			MOD	PPSM	MPPS for Group Case: SYS_NAME^YZ	R
420	CARD-2			MOD	IM	C-Store	R
450	NULL			--	--	Connectathon Test Completed!!!!	C

C-MOD IM DSS/OF CATH [Cardiology: CATH]

CATH Interactions Involving Modality, Image Manager, DSS/OF Trio

Ideally, there would be a hemodynamic (HD) and an angio (XA) modality available for these tests, however, that will most likely not be the case. If there is no HD modality available, execute with only a single modality for a cath lab. In the real world, this is analogous to the hemo system being non-IHE compliant.

This test primarily tests the DSS/OF ability to handle the complexities of CARD-1 transaction, but also tests the modalities ability to handle the various MPPS cases.

Modality performs different procedures (with different MPPS cases) and communicates those results with the Order Filler. Image Manager also participates. The MPPS Mgr can be combined with either the DSS/OF or the IM in these test cases, but only one MPPS Mgr should be active at one time and it should be recorded which MPPS Mgr was used.

We are testing the Modality worklist query. Show the observer that the modality can query the MWL using any combination of the following keys : Patient's Name (0010,0010), Patient ID (0010,0020), Accession Number (0008,0050), Requested Procedure ID (0040,1001). Demonstrate to the test observer that the DSS/OF is correctly creating and/or cancelling additional SPSs as necessary.

Note that the additional SPS defined below can be dynamically created by the DSS/OF if that is how the system is configured. The end result must be as shown, however. This test tests CATH Use Cases C1, C2 and C7.

Integration Profiles / Actors Concerned by the test :

<u>Integration Profile</u>	<u>Actor</u>	<u>Profile Option</u>	<u>Option</u>	<u>Meta Test</u>
CATH	IM	NONE	R	
CATH	MOD	NONE	R	
CATH	OF	NONE	R	

Test script :

Step	Trans.	Msg Type	Return	From Actor	To Actor	Step Description	Opt.
20	NULL			--	--	In the descriptions to follow, SYS_NAME is a combination of the Modality and DSS/Order Filer systems under test	R
30	NULL			--	--	For example: GEXAPHILIPSHDKODAKOF	R
40	NULL			--	--	You will need an DSS/OF for this test, but the ADT and OP may be the MESA RIS Mall.	R
50	NULL			--	--	You will need an Image Manager for this test; a simulator is not allowed	R
60	NULL			--	--	Use the ADT system to register 5 inpatients: SYS_NAME^AB, SYS_NAME^CD, SYS_NAME^EF,	R

						SYS_NAME^GH	
70	NULL			--	--	Use the Order Placer and place single orders for AB, CD, EF, GH	R
80	NULL			--	--	See the Table in the CATH Connectathon Test Plan for more explanation, if necessary.	R
90	NULL			--	--	Schedule 1 SPS XA and 1 SPS HD for SYS_NAME^AB	R
110	NULL			--	--	Schedule 1 SPS XA and 1 SPS HD for SYS_NAME^EF	R
120	NULL			--	--	Schedule 1 SPS XA and 1 SPS HD (different location) for SYS_NAME^GH	R
140	RAD-4			OF	IM	Scheduling messages for SYS_NAME^AB	R
160	RAD-4			OF	IM	Scheduling message for SYS_NAME^CD	R
180	RAD-4			OF	IM	Scheduling message for SYS_NAME^EF	R
200	RAD-4			OF	IM	Scheduling message for SYS_NAME^GH	R
220	NULL			--	--	Step 1: This begins the Simple 1 SPS:1MPPS Case for SYS_NAME^AB. The following transactions do not need to be in the exact order, just the sum total.	R
240	RAD-5			MOD	OF	Both modalities retrieve single SPS MWL for SYS_NAME^AB	R
260	CARD-1			MOD	PPSM	HD MPPS le in the exact order.	R
440	RAD-5			MOD	OF	Both modalities retrieve double SPS MWL for SYS_NAME^CD	R
870	CARD-1			MOD	PPSM	XA MPPS Intervention In Progress: SYS_NAME^EF	R
900	CARD-2			MOD	IM	XA C-Store	R
910	CARD-3			MOD	IM	HD Storage Commit	RO
911	CARD-3			MOD	IM	XA Storage Commit	R
920	RAD-7			MOD	PPSM	HD MPPS Intervention Complete	R
925	RAD-7			MOD	IM	XA MPPS Intervention Complete	R
930	NULL			--	--	Step 4: This begins the Broad Modality Query	R
940	RAD-5			MOD	OF	Both modalities retrieve single SPS	R

						MWL for SYS_NAME^GH	
950	CARD-1			MOD	PPSM	HD MPPS Diag In Progress: SYS_NAME^GH	R
960	CARD-1			MOD	PPSM	XA MPPS Diag In Progress: SYS_NAME^GH	R
970	NULL			--	--	Step 5: This begins the patient update portion of this test	R
980	NULL			OF	IM	Rename SYS_NAME^AB to SYS_NAME^AAABBB	R
990	NULL			OF	IM	Merge SYS_NAME^CD, SYS_NAME^EF to SYS_NAME^CD	R
999	NULL			--	--	Connectathon Test Completed!!!!	R

C-MOD IM DSS/OF ECHO [Cardiology: ECHO]

CATH Interactions Involving Modality, Image Manager, DSS/OF Trio

This test primarily tests the DSS/OF ability to handle the complexities of CARD-1 transaction, but also tests the modalities ability to handle the various MPPS cases.

Modality performs different procedures (with different MPPS cases) and communicates those results with the Order Filler. Image Manager also participates. The MPPS Mgr can be combined with either the DSS/OF or the IM in these test cases, but only one MPPS Mgr should be active at one time and it should be recorded which MPPS Mgr was used.

This test tests ECHO Use Cases E1, E2, E3, and E4.

Integration Profiles / Actors Concerned by the test :

<u>Integration Profile</u>	<u>Actor</u>	<u>Profile Option</u>	<u>Option</u>	<u>Meta Test</u>
ECHO	IM	NONE	R	
ECHO	MOD	NONE	R	
ECHO	OF	NONE	R	

Test script :

Step	Trans.	Msg Type	Return	From Actor	To Actor	Step Description	Opt.
20	NULL			--	--	In the descriptions to follow, SYS_NAME is a combination of the Modality and DSS/Order Filer systems under test	R
40	NULL			--	--	For example: GEECHO_KODAKOF	R
50	NULL			--	--	You will need an ADT and OP for this test. The MESA RIS Mail may be used for these actors, but not as the DSS/OF.	R
60	NULL			--	--	You will need an Image Manager for this test; a simulator is not allowed	R
80	NULL			--	--	Use the ADT system to register 5 inpatients: SYS_NAME^AA, SYS_NAME^BB, SYS_NAME^CC, SYS_NAME^DD	R
100	NULL			--	--	Use the Order Placer and place single orders for AA, BB, CC, DD	R
120	NULL			--	--	See the Table in the ECHO Connectathon Test Plan for more explanation, if necessary.	R
140	NULL			--	--	Schedule 1 SPS US for SYS_NAME^AA	R
160	NULL			--	--	Schedule 1 SPS US TEE and 1 SPS	R

						for US TTE for SYS_NAME^BB	
180	NULL			--	--	Schedule 1 SPS US (different location) for SYS_NAME^DD	R
200	RAD-4			OF	IM	Scheduling messages for SYS_NAME^AA	R
210	RAD-4			OF	IM	Scheduling messages for SYS_NAME^BB	R
220	RAD-4			OF	IM	Scheduling messages for SYS_NAME^DD	R
240	NULL			--	--	Step 1: This begins the Simple 1 SPS:1MPPS Case for SYS_NAME^AA. The following transactions do not need to be in the exact order, just the sum total.	R
260	RAD-5			MOD	OF	Echo modality retrieves single SPS MWL for SYS_NAME^AA	R
280	CARD-1			KER_SVR	PPSM	US MPPS In Progress: SYS_NAME^AA	R
300	CARD-2			MOD	IM	US C-Store	R
320	CARD-3			MOD	IM	US Storage Commit	R
500	CARD-2			MOD	IM	US C-Store	R
540	RAD-7			MOD	PPSM	US MPPS Intervention Complete for TTE	R
560	NULL			--	--	Step 3: This begins the Ad Hoc (add on, not scheduled) Case for patient SYS_NAME^CC. Start at the Modality this time, not with an order. Verify Modality Study Instance UID used.	R
580	CARD-1			MOD	PPSM	US MPPS Diag In Progress: SYS_NAME^CC	R
600	CARD-2			MOD	IM	US C-Store	R
620	CARD-3			MOD	IM	US Storage Commit	R
640	RAD-7			MOD	PPSM	US MPPS Diag Complete	R
660	NULL			--	--	Verify that the DSS/OF used the Study Instance UID of the modality.	R
680	NULL			--	--	Step 4: This begins the Broad Modality Query. This patient should be schedule for some other Echo modality/system or room. Force a broad query by the modality.	R
700	RAD-5			MOD	OF	Echo Modality retrieves single SPS	R

						MWL for SYS_NAME^DD	
720	CARD-1			MOD	PPSM	US MPPS Diag In Progress: SYS_NAME^DD	R
740	NULL			--	--	Step 5: This begins the patient update portion of this test which tests the IM	R
760	NULL			--	--	Rename SYS_NAME^AA to SYS_NAME^Aaaltered at ADT	R
780	NULL			--	--	Verify patient name change at IM.	R
800	NULL			--	--	Merge SYS_NAME^BB, SYS_NAME^CC to SYS_NAME^BB at ADT	R
820	NULL			--	--	Verify patients merged at IM.	R
840	NULL			--	--	Connectathon Test Completed!!!!	R

C-MOD IM ECHO [Cardiology: ECHO]

ECHO Profile: Modality/IM Communication

Modality performs different procedures (with different MPPS cases) and communicates those results with the Image Manager. This test assumes that the MPPS Mgr is associated with the IM. The IM/MMPS Mgr does not need to forward the MPPS messages to the OF for this test.

If you find a DSS/OF that wishes to participate, run the MOD_IM_DSS/OF test instead of this test. Note that this test case does NOT test the complexities of the DSS/OF for the CARD-1 transaction. See the MOD_IM_DSS/OF tests.

Integration Profiles / Actors Concerned by the test :

<u>Integration Profile</u>	<u>Actor</u>	<u>Profile Option</u>	<u>Option</u>	<u>Meta Test</u>
ECHO	IM	NONE	R	
ECHO	MOD	NONE	R	

Test script :

Step	Trans.	Msg Type	Return	From Actor	To Actor	Step Description	Opt.
20	NULL			--	--	In the descriptions to follow, SYS_NAME is a combination of the Modality and Image Mgr systems under test	R
30	NULL			--	--	NULL For example: GEECHO_SIEMENSIM	R
40	NULL			--	--	Find an ADT, Order Placer, DSS trio to help with this test. This may be the MESA RIS Mall.	R
50	NULL			--	--	Use the ADT system to register 5 inpatients: SYS_NAME^FF, SYS_NAME^GG, SYS_NAME^HH, SYS_NAME^JJ, SYS_NAME^KK	R
60	NULL			--	--	Use the Order Placer and place single orders for GG, HH, JJ	R
70	NULL			--	--	Use the Order Placer and place two orders for KK that can be used for the Group Case	R
80	NULL			--	--	Find a OF/MWL provider and schedule 1 SPS or SYS_NAME^GG	R
90	NULL			--	--	Find a OF/MWL provider and schedule 1 SPS for SYS_NAME^HH	R
100	NULL			--	--	Find a OF/MWL provider and schedule 1 SPS for SYS_NAME^JJ	R
110	NULL			--	--	Find a OF/MWL provider and schedule 2 SPSs for SYS_NAME^KK	R

120	NULL			--	--	This begins the PPS Unscheduled Case for patient SYS_NAME^FF; perform one (unscheduled) PPS for SYS_NAME^FF	R
130	CARD-1			MOD	PPSM	Echo MPPS for Unscheduled Case: SYS_NAME^FF	R
140	CARD-2			MOD	PPSM	C-Store acquisition data (images or reports)	R
150	CARD-2			MOD	IM	Storage Commitment	R
160	NULL			--	--	This begins the PPS Simple Case for patient SYS_NAME^GG; perform one (scheduled) PPS for SYS_NAME^GG	R
170	RAD-5			MOD	OF	MWL for SYS_NAME^GG	R
180	CARD-1			MOD	PPSM	MPPS for Simple Case: SYS_NAME^GG	R
190	CARD-2			MOD	IM	C-Store	R
200	CARD-3			MOD	IM	Storage Commitment	R
210	NULL			--	--	This begins the PPS Abandoned Case for patient SYS_NAME^HH; select this SPS and abandon the step. Abandon by sending an MPPS In Progress, then MPPS Discontinued.	R
220	CARD-1			MOD	PPSM	MPPS for Abandoned Case (In Progress then Cancel): SYS_NAME^HH	R
230	NULL			--	--	This begins the PPS Append Case for patient SYS_NAME^JJ; perform one (scheduled) PPS for SYS_NAME^JJ, and then append the same type of PPS	R
240	CARD-1			MOD	PPSM	MPPS for first step (In Progress then Complete) of Append Case: SYS_NAME^JJ	R
250	CARD-2			MOD	IM	C-Store	R
260	CARD-1			MOD	PPSM	MPPS for Append PPS (In Progress then Complete): SYS_NAME^JJ	R
270	CARD-2			MOD	IM	C-Store	R
280	NULL			--	--	This begins the PPS Group Case for patient SYS_NAME^KK; modality groups two SPS in to one PPS. Note that this step is optional for ECHO, only complete if the modality is	C

						capable of grouping procedures.	
290	RAD-5			MOD	OF	MWL for SYS_NAME^KK	R
300	CARD-1			MOD	PPSM	MPPS for Group Case (single In Progress/Complete for both SPS): C SYS_NAME^KK	
310	CARD-2			MOD	IM	C-Store	RO
320	NULL			--	--	Connectathon Test Completed!	R

C-OF Cancel [Cardiology: CATH]

CATH Use Case C8: Cancel at DSS/OF

The purpose of this test is to verify that an OF can cancel a procedure before it begins. This test only needs to be executed once per OF. This test tests CATH Use Case C8.

Integration Profiles / Actors Concerned by the test :

<u>Integration Profile</u>	<u>Actor</u>	<u>Profile Option</u>	<u>Option</u>	<u>Meta Test</u>
CATH	OF	NONE	R	
ECHO	OF	NONE	R	

Test script :

Step	Trans.	Msg Type	Return	From Actor	To Actor	Step Description	Opt.
20	RAD-2			OP	OF	or RAD-3 Create an order for a procedure.	R
30	NULL			--	--	Create a Requested Procedure/SPS at the OF	R
40	RAD-4			OF	IM	Send order notification to IM from OF. (Procedure Scheduled)	R
50	NULL			--	--	Cancel the Order at the OF	R
60	RAD-3			OF	OP	Cancel the order from the OF to the OP	R
70	RAD-13			OF	IM	Cancel the order from the OF to the IM. (Procedure Updated)	R
80	NULL			--	--	Note whether the cancelled procedure is noted at the Order Filler (preferred but not required).	O

ECG Adv Display [Cardiology: ECG]

Display Retrieves Multiple ECG Reports for display from Information Source

The ECG_Adv_Display test is used by the Display Actor to retrieve multiple reports for display from an Information Source. The Project Manager will observe the Display actor requesting each of the various types listed in the IHE Cardiology Technical Framework for this Integration Profile. The test will use one or more patients known to the Information Source, possibly those loaded during ECG_LoadInformationSource.

If the Information Source supports other Summary Request types, the RID test cases should be executed as well.

This test case tests the CARD-5 and CARD-6 Transactions.

For evaluation, the Project Manager will observe the query and result for multiple reports. Verify visually that the ECG page looks reasonable and contains the following information:

Patient Name, Patient ID, local date and time of ECG recording, voltage and timescale, lead labels, and frequency information, interpretation information and measurements . The ECG must state either CONFIRMED or UNCONFIRMED.

Integration Profiles / Actors Concerned by the test :

Integration Profile	Actor	Profile Option	Option	Meta Test
-------------------------------------	-----------------------	--------------------------------	------------------------	---------------------------

Test script :

Step	Trans.	Msg Type	Return	From Actor	To Actor	Step Description	Opt.
20	NULL			INFO_SRC	--	Determine one or more patients in Information Source with at least two reports.	R
40	CARD-5	HTTP Get		DISPLAY	INFO_SRC	Retrieve SUMMARY-CARDIOLOGY-ECG	R
60	CARD-6	HTTP Get		DISPLAY	INFO_SRC	Select and Retrieve two documents for the same patient	R
80	NULL			DISPLAY	--	Verify display of two ECGs	R

ECG LoadInformationSource [Cardiology: ECG]

ECG Information Sources Load Initial Data

The ECG_LoadInformationSource test is an initialization step for the Information Source to load initial data into its system. The Information Source can use the Patient Names and Patient IDs listed in this test or may choose other names and identifiers. The goal of the test is to have a known list of patients that are readily available for testing with Display actors.

There is no formal evaluation required. Tell the Project Manager you have completed this test when your system is properly loaded. You may find that you have to load new data later; do not report updated information.

The suggested patient list is:

- Peter Peterson u2013 234561
- Tommy Thompson u2013 234562
- Johnny Johnson - 234563

Integration Profiles / Actors Concerned by the test :

Integration Profile	Actor	Profile Option	Option	Meta Test
ECG	INFO_SRC	NONE	R	

Test script :

Step	Trans.	Msg Type	Return	From Actor	To Actor	Step Description	Opt.
20	NULL			INFO_SRC	--	Load initial data for patients listed or other patients.	R

ECG Retrieve [Cardiology: ECG]

Display Retrieves ECG Reports from Information Source

The ECG_Retrieve test is used by the Display Actor to retrieve various types of reports from an Information Source. The Project Manager will observe the Display actor requesting each of the various types listed in the IHE Cardiology Technical Framework for this Integration Profile. The test will use one or more patients known to the Information Source, possibly those loaded during ECG_LoadInformationSource.

If the Information Source supports other Summary Request types, the RID test cases should be executed as well. This test case tests the CARD-6 Transaction.

For evaluation, the Project Manager will observe the query and result for each summary type. Verify visually that the ECG page looks reasonable and contains the following information:

Patient Name, Patient ID, local date and time of ECG recording, voltage and timescale, lead labels, and frequency information, interpretation information and measurements . The ECG must state either CONFIRMED or UNCONFIRMED.

Integration Profiles / Actors Concerned by the test :

<u>Integration Profile</u>	<u>Actor</u>	<u>Profile Option</u>	<u>Option</u>	<u>Meta Test</u>
ECG	DISPLAY	NONE	R	
ECG	INFO_SRC	NONE	R	

Test script :

Step	Trans.	Msg Type	Return	From Actor	To Actor	Step Description	Opt.
20	NULL			INFO_SRC	--	Determine one or more patients in Information Source with reports.	R
40	CARD-5	HTTP Get		DISPLAY	INFO_SRC	Retrieve SUMMARY-CARDIOLOGY-ECG	R
60	CARD-6	HTTP Get		DISPLAY	INFO_SRC	Retrieve a document known to be of type Vector PDF.	R
80	CARD-6	HTTP Get		DISPLAY	INFO_SRC	Retrieve a document known to be of type SVG, if supported.	O
100	CARD-6	HTTP Get		DISPLAY	INFO_SRC	If the Display supports this option, retrieve a document directly using an OID (ie., do not perform a query and select from a list)	O

ECG SummaryRequest HTML [Cardiology: ECG]

ECG Summary Request HTML

The ECG_SummaryRequest_HTML test is used by the Display Actor to request summary reports from an Information Source. The Project Manager will observe the Display actor requesting each of the Summary Request types listed in the IHE Cardiology Technical Framework for this Integration Profile. The test will use one or more patients known to the Information Source, possibly those loaded during ECG_LoadInformationSource.

If the Information Source supports other Summary Request types, the RID test cases should be executed as well.

For evaluation, the Project Manager will observe the query and result for each summary type. Verify list of ECGs is returned via HTML page. Verify visually that the ECG page looks reasonable and contains the following information:

Patient Name, Patient ID, local date and time of ECG recording, voltage and timescale, lead labels, and frequency information, interpretation information and measurements . The ECG must state either CONFIRMED or UNCONFIRMED.

Integration Profiles / Actors Concerned by the test :

<u>Integration Profile</u>	<u>Actor</u>	<u>Profile Option</u>	<u>Option</u>	<u>Meta Test</u>
ECG	DISPLAY	NONE	R	
ECG	INFO_SRC	NONE	R	

Test script :

Step	Trans.	Msg Type	Return	From Actor	To Actor	Step Description	Opt.
20	NULL			INFO_SRC	--	Determine one or more patients in Information Source with reports.	R
40	ITI-11	HTTP Get		DISPLAY	INFO_SRC	Retrieve SUMMARY	R
60	ITI-11	HTTP Get		DISPLAY	INFO_SRC	Retrieve SUMMARY-CARDIOLOGY	R
80	NULL			DISPLAY	--	Select from list and display.	R

ECG SummaryRequest XML [Cardiology: ECG]

ECG Summary Request HTML

The ECG_SummaryRequest_XML test is used by the Display Actor to request summary reports from an Information Source. The Project Manager will observe the Display actor requesting each of the Summary Request types listed in the IHE Cardiology Technical Framework for this Integration Profile. The test will use one or more patients known to the Information Source, possibly those loaded during ECG_LoadInformationSource.

If the Information Source supports other Summary Request types, the RID test cases should be executed as well. This test case tests the CARD-5 Transaction- XML list responses.

For evaluation, the Project Manager will observe the query and result for each summary type. Verify list of ECGs is returned via XML list. Verify visually that the ECG page looks reasonable and contains the following information:

Patient Name, Patient ID, local date and time of ECG recording, voltage and timescale, lead labels, and frequency information, interpretation information and measurements . The ECG must state either CONFIRMED or UNCONFIRMED.

Integration Profiles / Actors Concerned by the test :

Integration Profile	Actor	Profile Option	Option	Meta Test
-------------------------------------	-----------------------	--------------------------------	------------------------	---------------------------

Test script :

Step	Trans.	Msg Type	Return	From Actor	To Actor	Step Description	Opt.
20	NULL			INFO_SRC	--	Determine one or more patients in Information Source with reports.	R
40	CARD-5	HTTP Get		DISPLAY	INFO_SRC	Retrieve SUMMARY-CARDIOLOGY-ECG	R
60	NULL			DISPLAY	--	Select from list and display.	R